

Vacuum Regulator High-Spatz -90 DIN Probe



Manufacturer: Greggersen

Product 900897

Number:

Description

Vacuum Regulator High-Spatz-90, Plug-in Unit DIN

Medical vacuum regulator with vacuum gauge,
fine dosing and quick-action shut-off valve.

Technical data:

Gas type: Vacuum / VAC

Inlet: Probe to DIN 13260 Part 2

Control range: 0 to -85 kPa (0 to -0.85 bar)

Outlet: 9/16-18 UNF with vacuum tubing nipple