


Forano, Base Part, DVE, DIN, O2



Manufacturer: Greggersen

Product 903230

Number:

Description

Terminal Unit Forano Base Part DIN/AGA/NF

Gas-specific base part with hose barb
for installation in ceiling-mounted supply systems.

Base part with mounting bracket, integrated non-return and maintenance valve.

Scope of delivery:

- gas-specific base part with hose barb
- mounting bracket

Technical data:

Gas type: Oxygen / O2

Installation type: DVE

Inlet: Hose barb 7 mm, rear

Properties

Gas Type O2-Oxygen