

Compressed-Air Vacuum Regulator Specht -90 Rail

NIST



Manufacturer: Greggersen

Product 902625

Number:

Description

Compressed-Air Vacuum Regulator Specht-90, Rail-Mounted Unit NIST

Medical compressed-air vacuum regulator for vacuum generation, with vacuum gauge, fine-adjustment and quick-shut-off valve, silencer and filter holder.

For mounting on an equipment rail (25x10 mm).

The gas inlet is integrated into the rail clamp (NIST).

Accessories required for operation:

Low-Pressure Hose AIR with NIST outlet

(not included in scope of delivery)

Technical Data:

Gas Type: Compressed Air / AIR

Inlet: NIST housing according to DIN EN ISO 18082

Control Range: 0 to -85 kPa (0 to -0.85 bar)

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	AIR Compressed Air
-----------------	--------------------
