


Falke without Flow O2 BS BULL +1BS



Manufacturer: Greggersen

Product 910344

Number:

Description

Falke without Flow with 1 quick coupling,
Compact pressure regulator with cylinder connection
according to BS 341-3:2002, contents gauge, with one
quick coupling at a pressure of 450 kPa.

Gas type: Oxygen / O2

Connection: High-pressure cylinder connection according to BS 341-3:2002

Connection no.: 3

Outlet: Quick coupling type BS 5682:1998, top

Output: 450kPa

Properties

Gas Type	O2-Oxygen
-----------------	-----------