

Falke w. Flow Stand. AIR DIN 9/16



Manufacturer: Greggersen

Product 910360

Number:

Description

AIR Falke with Flow, Standard,

Compact pressure regulator with contents gauge;
incrementally adjustable flow rate.

Gas type: Compressed air / AIR

Connection: Cylinder connection for gas cylinders according to DIN 477, No. 13

Outlet: 9/16-18 UNF (tubing nipple), bottom

Flow rate: 0...15 l/min

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type AIR Compressed Air