

Habicht Rail, O2, NIST, Kolibri



Manufacturer: Greggersen

Product 902676

Number:

Description

Emergency Supply Unit

Suction and insufflation unit, for oxygen

RAIL-MOUNTED UNIT NIST

Combi-unit for attachment to an equipment rail;
gas supply via a pressure regulator or
a terminal unit; ejector regulator with
vacuum gauge, fine control valve and silencer
with filter holder; Kolibri compact flowmeter with
10 flow settings

Gas type: Oxygen / O₂

Connection: NIST according to DIN EN ISO 18082

Flowmeter outlet: 9/16-18 UNF with tubing nipple

Ejector outlet: Vacuum tubing nipple

consisting of:

- Rail clamp with pressure inlet NIST
- Distribution block, spacer
- Kolibri standard 0...15 l/min
- Ejector regulator, 0 to -0.9 bar