

Kolibri 2 Probe MEDAP AIR 0/8



Manufacturer: Greggersen

Product 906455

Number:

Description

KOLIBRI 2, AIR, OFF / 8 l/min

- MEDAP PLUG-IN UNIT -

Flow Settings: 0 / 8 l/min

Compact Flowmeter KOLIBRI 2
with fixed orifices with selector device

Medical compact flowmeter for coupling
into a terminal unit, with handwheel - stepwise
adjustable flow rate, outlet 9/16-18 UNF thread

Ergonomically optimised adjustment device and readability
of the selected flow setting from various viewing directions
and angles.

Frontal operating concept -> higher safety during
setting of the correct flow rate

Gas type: med. compressed air / AIR
Connection: MEDAP probe
Outlet: 9/16-18 UNF
Flow rate: 8 l/min at inlet pressure 450 kPa (4.5 bar)

Accessories required for operation:
900619 - Hose connection 9/16-18 UNF

Properties

Gas Type	AIR Compressed Air
Inlet Standard	MEDAP
Variant	Air8