

Vacuum Regulator Spatz -16 Rail NIST



Manufacturer: Greggersen

Product 902624

Number:

Description

Vacuum Regulator SPATZ-16, VAC, Rail-Mounted Unit NIST

Medical vacuum regulator, with vacuum gauge,
fine dosing and quick-shut-off valve,
for mounting on an equipment rail (25x10mm)

Accessories required for operation:

Low-Pressure Hose VAC with NIST outlet
(not included in scope of delivery)

Technical data:

Gas type: Vacuum

Connection: NIST housing according to DIN EN ISO 18082

Outlet: Vacuum tubing nipple

Output: adjustable between 0 and -16 kPa (0 to -160 mbar)