

Vacuum Reg. Skua 2.0 -30,AIR,RAIL,NIST



Manufacturer: Greggersen

Product 904603

Number:

Description

Compressed-Air Vacuum Regulator Skua -30, Rail-Mounted Unit NIST

Medical compressed-air vacuum regulator for dosing vacuum.

The vacuum is generated by means of compressed air using the Venturi principle.

With 360° rotatable vacuum gauge; large, ergonomic adjustment handwheel; rotatable exhaust sleeve.

The gas inlet is integrated into the rail clamp (NIST).

The adjustment mechanism enables stepless, precise, and intuitive vacuum control. The rotatable vacuum gauge allows the user to read the set value from an optimal viewing angle.

The exhaust sleeve also allows the air flow to be aligned individually.

The regulator unit complies with the requirements of DIN EN ISO 10079-3:2014 (high vacuum /high flow).

Material:

Housing: aluminium, anodised

Handwheel: plastic

Rail clamp: aluminium, anodised

Technical Data:

Gas type: med. compressed air / AIR

Inlet: NIST housing in accordance with DIN EN ISO 18082

Control range: 0 to -30 kPa $\pm$ 10%

Flow rate: >20 l Freeflow/min

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	AIR Compressed Air
-----------------	--------------------
