

Kolibri, DIN Probe, O2 Standard+ long



Manufacturer: Greggersen

Product 901481

Number:

Description

Kolibri, O₂, 0...30 l/min, DIN

- PLUG-IN UNIT, long -

Medical compact flowmeter in accordance with DIN EN ISO 15002,
for coupling into a terminal unit, with handwheel
step-adjustable flow rate, outlet 9/16-18 UNF
and tubing nipple

Gas type: Oxygen / O₂

Connection: Probe according to DIN 13260 Part 2

Total length: Probe front part to centre of flow outlet 108 mm

Outlet: 9/16 - 18 UNF / tubing nipple

Flow rate: 0...30 l/min

Flow settings: 0/ 0.5/ 1/ 2/ 4/ 6/ 8/ 10/ 12/ 15/ 30+l/min

Properties

Gas Type	O ₂ -Oxygen
-----------------	------------------------
