

Ventus evo 6-way



Manufacturer:	Greggersen
----------------------	------------

Product	905206
----------------	--------

Number:	
----------------	--

Description

Ventus evo – 6-way Area Control Unit for medical gases and vacuum

VENTUS evo is an Area Control Unit in accordance with DIN EN ISO 7396-1 and HTM 02-01

for reliable monitoring and control of medical gases and vacuum.

Designed to accommodate six valve inserts (available separately).

Ensures safe transmission of the respective media.

Integrated pressure monitoring is carried out via sensors; in addition, there is the option for shut-off in case of emergency or installation, as well as emergency supply via a NIST connection.

Monitoring of operational statuses and the output of alarm and warning messages are performed via the integrated area monitor. An on-board relay for common fault indication is integrated as standard.

Technical Data

Mains connection: 100–240 VAC, 50/60 Hz

Design: Flush-mounted / cavity wall mounting

Plaster compensation: max. 25 mm

Dimensions Front Part (W × H × D): 630 × 517 × 17 mm

Dimensions Base Part (W × H × D): 568 × 470 × 68.5 mm

Medical Device: Class IIb

Configuration

Gas 1: ### please specify! ###

Gas 2: ### please specify! ###

Gas 3: ### please specify! ###

Gas 4: ### please specify! ###

Gas 5: ### please specify! ###

Gas 6: ### please specify! ###

Language: German

Pressure unit: bar

Sensor inputs: -1 to +6 bar

Alarm limits O₂, AIR, N₂O, CO₂, N₂: 500 kPa ± 20 %

Alarm limit MA4: 400 kPa ± 20 %

Alarm limit VAC: -40 kPa

Alarm limit AIR-800: 800 kPa ± 20 %

Alarm limit SA7: 700 kPa ± 20 %

Door alarm: yes / no

Relay module: no / 2

Modbus module: yes / no

Network: activated / not activated