

Kolibri, DIN Probe, AIR 10 Flows



Manufacturer: Greggersen

Product 901407

Number:

Description

Kolibri, Air, 10 Flows

-Plug-in Unit-

Medical Compact Flowmeter in accordance with
DIN EN ISO 15002, for coupling
into a terminal unit, with handwheel for
incrementally adjustable flow rate,
outlet 9/16-18 UNF and tubing nipple

Gas type: Compressed Air

Connection: Probe according to DIN 13260 Part 2

Total length: 88 mm (Probe Front Part to centre of Tubing Nipple outlet)

Outlet: 9/16-18 UNF / Tubing Nipple

Flow rate: 0... l/min

Flow settings: (Please specify flow settings!)

Properties

Gas Type	AIR Compressed Air
----------	--------------------
