


Forano Front Part, DIN, N2-800, Colour-Neutral



Manufacturer: Greggersen

Product 905116

Number:

Description

Terminal Unit FORANO - Front Part DIN

Colour-Neutral

Gas type: Nitrogen 800 / N2-800

Labelling: nitrogen 800 / N2-800

Gas-specific socket assembly for mounting
on FORANO base part.

Socket assembly according to DIN 13260 Part 2
with parking and operating position.

Release via release ring.

Colour coding:

Colour-Neutral (black)

Release ring labelling: printed

Language: English

Cartridge incl. sealing elements

Technical Data:

Gas type: Nitrogen 800 / N2-800

Labelling: nitrogen 800 / N2-800

Colour coding: Colour-Neutral (black)

Front diameter: Ø 42mm

Scope of delivery:

- gas-specific socket assembly
- release ring incl. printed labelling
- cartridge with sealing elements

Properties

Gas Type	N2800-Nitrogen 8bar
-----------------	---------------------
