


Forano, Front Part, DIN, VAC, ISO



Manufacturer: Greggersen

Product 905102

Number:

Description

Terminal Unit FORANO - Front Part DIN

ISO 32 - Version III from 08/2019

Gas type: Vacuum / VAC

Inscription: vacuum / VAC

Gas-specific socket assembly for mounting
on base part FORANO.

Socket assembly according to DIN 13260 Part 2
with parking and operating position.

Unlocking via release ring.

Colour coding:

according to DIN EN ISO 9170-1 (ISO32 - yellow)

Release ring inscription: printed

Language: English

Cartridge incl. sealing elements

Technical Data:

Gas type: Vacuum / VAC

Inscription: vacuum / VAC

Colour coding: according to ISO 9170-1 (yellow)

Front diameter: Ø 42mm

Scope of Delivery:

- gas-specific socket assembly
- release ring incl. printed inscription
- cartridge with sealing elements

Properties

Gas Type	VAC-Vacuum
-----------------	------------
