


Falke with Standard Flow O2 DIN Tubing Nipple +1DIN



Manufacturer: Greggersen

Product 910204

Number:

Description

Falke with Flow, standard with 1 quick-release coupling,
Compact Pressure Regulator with contents gauge;
stepwise adjustable flow rate,
DIN quick-release coupling at a pressure of 450 kPa.

Gas type: Oxygen / O₂

Connection: Cylinder Connection for gas cylinders according to DIN 477, No. 9

Outlet: Tubing Nipple / 1 quick-release coupling according to DIN 13260 Part 2

Output: 0...15 l/min Tubing Nipple - 120 l/min quick-release coupling

Flow Settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type

O₂-Oxygen