

Conport V2A Pirol 2.0 -90, VAC, NIST



Manufacturer: Greggersen

Product 904611

Number:

Description

Portable ConPort Suction Unit with Pirol -90 Vacuum Regulator

Medical, portable suction unit - vacuum-driven -
including stainless steel carrying unit to accommodate 2
MEDELA 1.5 L or Serres 1.0 L suction jars.

The carrying frame features a sturdy handle and a
bracket for a 25x10 mm equipment rail.

The vacuum control unit is located at the front,
allowing the operator quick and direct access.

The gas inlet as a NIST coupling is located at the rear
of the unit, keeping the supply hose outside
the working area.

Scope of delivery:

- Pirol-90 Vacuum Regulator with fine dosing valve,
rotatable vacuum gauge
- Carrying frame, V2A stainless steel

- 2 holders for secretion jars

Technical Data:

Gas type: Vacuum - VAC

Connection: NIST according to DIN EN ISO 18082

Outlet: 9/16-18 UNF with tubing nipple

Control range: 0 to -90 kPa $\pm$ 10%

Weight: approx. 3 kg (without jars)

Properties

Gas Type

VAC-Vacuum
