


Hose, CO2, 5m, DIN - DIN Coupling Colour-Neutral



Manufacturer: Greggersen

Product 902538

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe with reverse flow limiter,
hose and coupling according to DIN

Specification:

Gas type: Carbon dioxide / CO2

Inlet: Probe according to DIN 13260 Part 2

Hose length: 5.0 m

Colour coding: Colour-Neutral

Outlet: Coupling according to DIN 13260 Part 2

Properties

Gas Type CO2-Carbon Dioxide