

Nebuliser Plug-in Valve with Socket Assembly AIR



Manufacturer: Greggersen

Product 900774

Number:

Description

Fitis Nebuliser Plug-in Valve with Socket AIR

Plug-in Unit

Nebuliser plug-in valve, preset with socket (top) and
Combi-Probe, for connection to a terminal unit.

Technical Data:

Gas Type: Compressed Air

Connection: Probe to DIN 13260 Part 2

Outlet: 9/16-18 UNF

Flow Rate: 5 l/min $\pm$ 20%

Properties

Gas Type AIR Compressed Air