


Ventus 22mm 6-way Front Part+Base Part



Manufacturer: Greggersen

Product 900882

Number:

Description

Floor Shut-off Valve Box VENTUS 22mm

6-way, Front Part + Base Part

Gas type: *** Please specify ***

Language of Instructions for Use: German, English, Dutch, Danish, Hungarian (please delete where not applicable!)

Design complies with DIN EN ISO 7396-1 and HTM 02-01.

Front part assembly consisting of:

- Sheet steel door, powder-coated
- Side opening via hinges
- 2-way alarm unit
- Magnetic foil to cover the warning
- Plaster compensation up to 25mm
- Power supply unit 230 V AC

Base part assembly consisting of:

- Built-in housing for flush mounting

The alarm unit monitors:

- Mains pressure too high
- Mains pressure too low

Dimensions

Inlet & Outlet: 22 mm copper pipe

Front Part (WxHxD): 630 x 530 x 9 mm

Base Part (WxHxD): 570 x 470 x 77 mm

Warning alarm limits:

O₂, AIR, N₂O, CO₂: 500 kPa ± 20%

VAC: -40 kPa

AIR-800: 800 kPa ± 20%