


Kolibri, AGA Probe, O2 Standard



Manufacturer: Greggersen

Product 901420

Number:

Description

Kolibri, O2, Standard

PLUG-IN UNIT

Medical compact flowmeter according to DIN EN ISO 15002,
for coupling into a terminal unit, with handwheel
incrementally adjustable flow rate, outlet 9/16-18 UNF
and tubing nipple

Gas: Oxygen

Connection: Probe Type SS 87 524 30 (AGA)

Outlet: 9/16-18 UNF / Tubing Nipple

Flow Rate: 0...15 l/min

Flow Settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type	O2-Oxygen
----------	-----------
