

Nebuliser Plug-in Valve with Clamp O2



Manufacturer: Greggersen

Product 900776

Number:

Description

Fitis Nebuliser Plug-in Valve with Clamp O2

Plug-in Unit

Nebuliser plug-in valve, preset with side clamp;
for coupling into a terminal unit

Technical Data:

Gas type: Oxygen

Connection: Probe to DIN 13260 Part 2

Outlet: 9/16-18 UNF/Tubing Nipple

Flow Rate: 5 l/min $\pm$ 20%

Properties

Gas Type	O2-Oxygen
-----------------	-----------