

Portable Emergency Supply Habicht Kolibri



Manufacturer: Greggersen

Product 902697

Number:

Description

Emergency Supply Unit,
portable, without oxygen cylinder
- with KOLIBRI Compact Flowmeter -

Emergency supply unit, oxygen flowmeter
and ejector for generating a vacuum
consisting of:

Stainless steel carrying frame with stand, to accommodate
a 2 l or 3 l O₂ cylinder, suitable
for mounting on equipment rails or bed,
central supply block with NIST inlet,
consisting of:
KOLIBRI Compact Flowmeter 0...15 l/min with flow settings
for flow dosing and ejector with vacuum gauge
and fine dosing valve;

Gas type: Oxygen / O₂

Connection: NIST according to DIN EN ISO 18082

Outlet: Tubing Nipple / Vacuum Nipple

Output: FLOW : 0...15 l/min

VAC: adjustable between

0 and -85 kPa (0 to -0.85 bar)