


AGSS Forano VAC DVE Type 1L



Manufacturer: Greggersen

Product 902078

Number:

Description

AGSS FORANO Type 1L Vacuum
for ceiling pendants

The AGSS FORANO Type 1L Vacuum assembly serves as a socket unit to accommodate an AGSS Type 1L probe and enables integration into an AGSS system in accordance with DIN EN ISO 9170-2.

Equipped with a robust connection nipple (Ø 15 mm), the module can be connected directly to a vacuum drive unit in accordance with EN ISO 7396-2.

Technical data:

Flow rate: ≥ 90 l/min

Inlet: Tubing Nipple Ø 15 mm

Gas: Vacuum

Properties

Gas Type	AGSS-UTL
----------	----------
