

# Compressed-Air Vacuum Regulator Specht -10 Rail NIST



**Manufacturer:** Greggersen

**Product** 902672

**Number:**

## Description

Compressed-Air Vacuum Regulator Specht-10, Rail-Mounted Unit NIST

Medical compressed-air vacuum regulator for vacuum generation,  
with vacuum gauge, fine adjustment and quick-acting valve,  
silencer and filter holder.

For mounting to an equipment rail (25x10 mm).

The gas inlet is integrated in the rail clamp (NIST).

Accessories required for operation:

Low-Pressure Hose AIR with NIST outlet

(not included in the scope of delivery)

Technical data:

Gas type: Compressed air / AIR

Inlet: NIST housing according to DIN EN ISO 18082

Control range: 0 to -10 kPa (0 to -100 mbar)

Outlet: 9/16-18 UNF with vacuum tubing nipple

## Properties

|                 |                    |
|-----------------|--------------------|
| <b>Gas Type</b> | AIR Compressed Air |
|-----------------|--------------------|

---