


Forano, Front Part, DIN, O2, ISO



Manufacturer: Greggersen

Product 905100

Number:

Description

Terminal Unit FORANO - Front Part DIN

ISO 32 - Version III from 08/2019

Gas type: Oxygen / O₂

Labelling: oxygen / O₂

Gas-specific socket assembly for mounting
on base part FORANO.

Socket assembly according to DIN 13260 Part 2
with parked and operating position. Release via release ring.

Colour coding:

according to DIN EN ISO 9170-1 (ISO32 - white)

Release ring labelling: printed

Language: English

Cartridge incl. sealing elements

Technical data:

Gas type: Oxygen / O₂

Labelling: oxygen / O₂

Colour coding: according to ISO 9170-1 (white)

Front diameter: Ø 42mm

Scope of delivery:

- gas-specific socket assembly
- release ring incl. printed labelling
- cartridge with sealing elements

Properties

Gas Type

O₂-Oxygen
