


Forano, Front Part, AGA, AIR, ISO



Manufacturer: Greggersen

Product 905141

Number:

Description

Terminal Unit FORANO - Front Part AGA

Gas type: Compressed Air / AIR

Inscription: medical air / AIR

Gas-specific socket assembly for mounting
on base part FORANO. Design complies with
DIN EN ISO 9170-1;

Socket assembly according to SS 87 524 30

Release via release ring.

Colour coding according to ISO 32.

Release ring inscription: printed

Cartridge incl. sealing elements

Technical Data:

Gas type: Compressed Air / AIR

Inscription: medical air / AIR

Colour coding: ISO 32

Front diameter: Ø 42mm

Scope of Delivery:

- gas-specific socket assembly
- release ring incl. printed inscription
- cartridge with sealing elements

Properties

Gas Type	AIR Compressed Air
----------	--------------------
