


Falke with Flow Stand CO2 DIN 9/16 + 1DIN



Manufacturer: Greggersen

Product 910354

Number:

Description

Falke with Flow, Standard with 1 quick-release coupling
 Compact pressure regulator with cylinder contents gauge;
 incrementally adjustable flow rate,
 DIN quick-release coupling at a pressure of 450 kPa.

Gas type: Carbon dioxide / CO2

Connection: Cylinder connection for gas cylinders according to DIN 477, No. 6

Outlet: 9/16-18 UNF (tubing nipple) / 1 quick-release coupling according to DIN 13260-2

Output: 0...15 l/min 9/16-18 UNF; 120 l/min quick-release coupling

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type	CO2-Carbon Dioxide
-----------------	--------------------