


Falke Flow Stand. O2 DIN Tubing Nipple +1DIN can



Manufacturer: Greggersen

Product 910212

Number:

Description

Falke with Flow, Standard

with 1 quick coupling

Compact pressure regulator with contents gauge;

incrementally adjustable flow rate,

DIN quick coupling at a pressure of 450 kPa.

Gas type: Oxygen / O₂

Connection: Cylinder connection for gas cylinders DIN 477, No. 9

Outlet: Tubing nipple, upwards

Outlet: Quick coupling DIN 13260-2, downwards

Flow rate: 0...15 l/min tubing nipple / 120 l/min quick coupling

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type

O2-Oxygen
