


AGSS Forano AIR DVE1 Type1L



Manufacturer: Greggersen

Product 902074

Number:

Description

AGSS FORANO Type 1L

for ceiling pendants with external ejector

according to DIN EN ISO 9170-2, with exhaust return,
visual indicator for operational control and with
mounting bracket and outlet pipe, flexible
hose connection to the ejector

Output: min. 50 l/min at 500 kPa (5 bar)

Gas type: Compressed Air / AIR

Inlet: Tubing Nipple Ø 7 mm

Outlet: Exhaust air Tubing Nipple Ø 15 mm

Drive gas Tubing Nipple Ø 7 mm

Properties

Gas Type

AIR Compressed Air
