


Pressure Regulator Unit Tank Supply



Manufacturer: Greggersen

Product 327800

Number:

Description

Pressure Regulator Unit for Tank Supply

Version with pressure sensor

Gas type: ***** (Please specify!)

Inlet pressure: max. 2,000 kPa (20 bar)

Flow rate: max. 150 m³/h;

Outlet pressure: 500 kPa (5 bar)

For reducing the pressure of a liquid gas tank;

suitable for use in medical supply systems

according to DIN EN ISO 7396-1. Max. inlet pressure 20 bar, outlet pressure 5 bar,

with shut-off valve, pressure regulator, pressure gauge, safety valve.

Completely mounted on mounting plate.

Inlet & outlet: 22mm solder connection

SV outlet: 15mm solder connection