

Flowmeter, NIST, AIR, 0-15l/min



Manufacturer: Greggersen

Product 902613

Number:

Description

Flowmeter, AIR, 0...15 l/min

RAIL-MOUNTED UNIT NIST

Medical flowmeter, pressure-compensated

in accordance with DIN EN ISO 15002, with fine adjustment valve

for dosing the flow; gas inlet

NIST integrated in rail clamp

Gas type: Compressed Air AIR

Inlet: NIST according to DIN EN ISO 18082

Outlet: 9/16-18 UNF

Flow Rate: 0...15 l/min

Properties

Gas Type

AIR Compressed Air