


Forano, Base Part, DVE, DIN, BS, VAC



Manufacturer: Greggersen

Product 903232

Number:

Description

Terminal Unit Forano Base Part DIN/AGA/NF/BS

Gas-specific base part with hose barb
for installation in ceiling supply systems.

Base part with mounting bracket, integrated non-return and maintenance valve.

Scope of delivery:

- Gas-specific base part with hose barb
- Mounting bracket

Technical data:

Gas type: Vacuum / VAC

Standard AGA: UTL (AGSS)

Installation type: DVE

Inlet: Hose barb 8.8 mm, rear

Properties

Gas Type

VAC-Vacuum
