


Falke without Flow easy 6l O2 DIN 9/16



Manufacturer: Greggersen

Product 910371

Number:

Description

Falke "easy" with flow,
Compact pressure regulator with contents gauge;
fixed flow rate of 6 l/min at a pressure of 450 kPa.

Gas type: Oxygen / O₂

Connection: Cylinder connection for gas cylinders to DIN 477, No. 9

Outlet: 9/16-18 UNF with tubing nipple at the bottom

Flow rate: 6 l/min

Properties

Gas Type O₂-Oxygen