

Kolibri 2 DIN Probe CO2 Standard



Manufacturer: Greggersen

Product 906419

Number:

Description

KOLIBRI 2, CO₂, 0...15 l/min

- Standard DIN PLUG-IN UNIT -

Flow Settings:

0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Compact Flowmeter KOLIBRI 2

with fixed orifices and selector device

Medical compact flowmeter for coupling
into a terminal unit, with handwheel – incrementally
adjustable flow rate, outlet 9/16-18 UNF thread

Ergonomically optimised adjustment device and readability
of the set flow setting possible from various viewing directions
and angles.

Frontal operating concept -> greater safety when
setting the correct flow rate

Gas type: Carbon dioxide / CO₂

Connection: Probe according to DIN 13260 Part 2

Outlet: 9/16-18 UNF

Flow rate: 0...15 l/min at inlet pressure 450 kPa (4.5 bar)

Accessories required for operation:

900619 - Hose connection 9/16-18 UNF

Properties

Gas Type	CO ₂ -Carbon Dioxide
Inlet Standard	DIN
Variant	Standard