

Double Kolibri, DIN, O2 Standard, ISO32



Manufacturer: Greggersen

Product 902163

Number:

Description

O2 - Double Kolibri, Standard

Medical compact flowmeter according to DIN EN ISO 15002,
for coupling into a terminal unit, with handwheel
incrementally adjustable flow rate, outlet 9/16-18 UNF
and tubing nipple

Handwheel in ISO32 - white for O2

PLUG-IN UNIT

Gas type: Oxygen

Connection: Probe according to DIN 13260 Part 2

Outlet: 9/16-18 UNF / tubing nipple

Flow rate: 2 x 0...15 l/min

Flow settings: 0/ 0.5/ 1/ 1.5/ 2/ 3/ 4/ 6/ 8/ 12/ 15 l/min

Properties

Gas Type	O2-Oxygen
----------	-----------
