

Habicht Rail, O2, NIST, Kolibri+Coupling



Manufacturer: Greggersen

Product 902677

Number:

Description

Emergency Supply Unit

Suction and insufflation unit, for oxygen

RAIL-MOUNTED UNIT NIST

Combi-unit for attachment to an equipment rail

Gas supply via a pressure regulator or
a gas terminal unit; ejector regulator with
vacuum gauge, fine control valve and silencer
with filter mount, Kolibri compact flowmeter with
10 flow settings and additional 5 bar outlet
as central gas supply coupling.

Gas type: Oxygen / O2

Connection: NIST according to DIN EN ISO 18082

consisting of:

- Rail clamp with NIST pressure inlet
- Distribution block, spacer
- Kolibri Standard 0...15 l/min
- Ejector regulator 0 to -0.9 bar
- Quick-connect coupling according to DIN 13260 Part 2