


Falke without Flow O2 NL BENE +1NF



Manufacturer: Greggersen

Product 500077

Number:

Description

Falke without Flow
with an NF quick-release coupling

Compact pressure regulator with high-pressure cylinder connection nut, contents gauge, quick-release coupling at a pressure of 450 kPa

Gas type: Oxygen

Connection: High-pressure cylinder connection

NF E29-650:1992, F-627

Outlet: NF quick-release coupling

Output: 450 kPa

Properties

Gas Type

O2-Oxygen