

Kolibri, NF Probe, O2 Standard



Manufacturer: Greggersen

Product 901460

Number:

Description

Kolibri, O2, Standard, NF

PLUG-IN UNIT

Medical compact flowmeter according to DIN EN ISO 15002,
for connecting to a terminal unit, with handwheel
gradually adjustable flow rate,
outlet 9/16-18 UNF and tubing nipple

Gas Type: Oxygen

Connection: Probe type NORME FRANCAISE

Outlet: 9/16-18 UNF / tubing nipple

Flow Rate: 0...15 l/min

Flow Settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type	O2-Oxygen
----------	-----------
