


Angled Probe CO2 Hose Barb 7mm, DIN, ISO



Manufacturer: Greggersen

Product 900639

Number:

Description

Angled Probe for Carbon Dioxide
with 7 mm Hose Barb and Non-Return Valve

Gas Type: Carbon Dioxide / CO₂

Connection: Probe according to DIN 13260 Part 2

Outlet: 7 mm Hose Barb

Colour Coding: according to ISO 9170-1 (grey)

Properties

Gas Type	CO ₂ -Carbon Dioxide
-----------------	---------------------------------