

Double Kolibri, DIN, AIR Standard



Manufacturer: Greggersen

Product 500248

Number:

Description

Double Kolibri - AIR

Standard/ PLUG-IN UNIT

Medical Compact Flowmeter according to DIN EN ISO 15002,
for coupling into a terminal unit, with handwheel
gradually adjustable flow rate, outlet 9/16-18 UNF
and tubing nipple

Gas type: Compressed Air - AIR

Connection: DIN Probe according to DIN 13260 Part 2

Outlet: 9/16-18 UNF / tubing nipple

Flow Rate: 2 x 0...15 l/min

Flow Settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type

AIR Compressed Air
