

Kolibri 2 Rail NIST CO2 Standard



Manufacturer: Greggersen

Product 906435

Number:

Description

KOLIBRI 2, CO₂, 0...15 l/min

- Standard RAIL-MOUNTED UNIT -

Flow Settings:

0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Compact Flowmeter KOLIBRI 2

with permanently attached orifices with selector mechanism

Medical compact flowmeter for mounting on an
equipment rail 25x10mm; gas inlet NIST integrated into
the rail clamp; with handwheel - step-adjustable
flow rate, outlet 9/16-18 UNF thread

Ergonomically optimised adjustment mechanism and readability
of the set flow setting from various viewing directions

and angles possible.

Front-facing operating concept -> increased safety when setting the correct flow rate

Gas type: Carbon dioxide / CO₂

Connection: NIST in accordance with DIN EN ISO 18082

Outlet: 9/16-18 UNF

Flow rate: 0...15 l/min at inlet pressure 450 kPa (4.5 bar)

Accessories required for operation:

900619 - Hose connection 9/16-18 UNF

Properties

Gas Type	CO ₂ -Carbon Dioxide
Inlet Standard	NIST
Variant	Standard