


Falke with FLOW Standard O2 DIN Tubing Nipple



Manufacturer: Greggersen

Product 910200

Number:

Description

Falke with Flow, Standard,
Compact pressure regulator with contents gauge;
step-adjustable flow rate.

Gas type: Oxygen / O₂

Connection: Cylinder connection for gas cylinders according to DIN 477, no. 9

Outlet: Tubing nipple below

Flow rate: 0...15 l/min

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type O₂-Oxygen