


Angled Probe CO2 Hose Barb 7mm, DIN,ISO V4A



Manufacturer: Greggersen

Product 505154

Number:

Description

Angled Probe for Carbon Dioxide
with 7 mm Hose Barb and Non-Return Valve
Probe Front Part made of 1.4404 (Stainless Steel)

Gas Type: Carbon Dioxide / CO2
Connection: Probe to DIN 13260 Part 2
Outlet: 7 mm Hose Barb
Colour Coding: to ISO 9170-1 (grey)

Properties

Gas Type	CO2-Carbon Dioxide
-----------------	--------------------