

Compressed-Air Vacuum Regulator Specht -10 DIN

Probe



Manufacturer: Greggersen

Product 901914

Number:

Description

Compressed-Air Vacuum Regulator Specht-10, Plug-in Unit DIN

Medical compressed-air vacuum regulator for vacuum generation, with vacuum gauge, fine adjustment and quick-action shut-off valve, silencer and filter holder.

Technical Data:

Gas Type: Compressed Air / AIR

Inlet: Probe according to DIN 13260 Part 2

Control Range: 0 to -10 kPa (0 to -100 mbar)

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type AIR Compressed Air