

Kolibri 2 Rail NIST AIR Standard+



Manufacturer: Greggersen

Product 906427

Number:

Description

KOLIBRI 2, AIR, 0...30 l/min

- Standard+ RAIL-MOUNTED UNIT -

Flow Settings:

0 / 2 / 3 / 4 / 6 / 8 / 12 / 15 / 20 / 25 / 30 l/min

Compact Flowmeter KOLIBRI 2

with fixed orifices and selection device

Medical compact flowmeter for mounting on an equipment rail 25x10mm; NIST gas inlet integrated into the rail clamp; with handwheel - incrementally adjustable flow rate, outlet 9/16-18 UNF thread

Ergonomically optimised adjustment device and legibility of the set flow setting possible from various viewing directions and angles.

Frontal operating concept -> greater safety when
setting the correct flow rate

Gas Type: Med. Compressed Air / AIR
Connection: NIST according to DIN EN ISO 18082
Outlet: 9/16-18 UNF
Flow Rate: 0...30 l/min at inlet pressure 450 kPa (4.5 bar)

Accessories required for operation:
900619 - Hose Connection 9/16-18 UNF

Properties

Gas Type	AIR Compressed Air
Inlet Standard	NIST
Variant	Standard+