

## Vacuum Reg. Pirol 2.0 -30,VAC,RAIL,NIST



**Manufacturer:** Greggersen

**Product** 904607

**Number:**

### Description

Vacuum Regulator Pirol -30, Rail-Mounted Unit NIST

Medical vacuum regulator for dosing vacuum.

With 360° rotatable gauge; large, ergonomic adjustment handwheel. The adjustment mechanism enables continuous, precise and intuitive vacuum regulation. The rotatable gauge allows the user to read the set value from an optimal viewing angle.

The regulator unit meets the requirements of DIN EN ISO 10079-3 (high vacuum / high flow).

Material:

Housing: aluminium, anodised

Handwheel: plastic

Rail clamp: aluminium, anodised

Technical data:

Gas type: Vacuum / VAC

Inlet: NIST housing according to DIN EN ISO 18082

Control range: 0 to -30 kPa  $\pm$  10%

Flow rate: > 20 l freeflow/min

Outlet: 9/16-18 UNF with vacuum tubing nipple

## Properties

| Gas Type |
|----------|
|----------|

|            |
|------------|
| VAC-Vacuum |
|------------|

---