


Angled Probe CO2 Hose Barb 7mm, DIN, Colour-Neutral



Manufacturer: Greggersen

Product 900608

Number:

Description

Angled Probe for Carbon Dioxide
with 7 mm hose barb and Non-Return Valve

Gas type: Carbon Dioxide / CO2

Connection: Probe to DIN 13260 Part 2

Outlet: 7 mm Hose Barb

Colour coding: neutral (black)

Properties

Gas Type	CO2-Carbon Dioxide
-----------------	--------------------