

Flowmeter, NIST, O2, 0-15l/min



Manufacturer: Greggersen

Product 902611

Number:

Description

Flowmeter, O2, 0...15 l/min

RAIL-MOUNTED UNIT NIST

Medical flowmeter, pressure-compensated

according to DIN EN ISO 15002, with fine regulating valve

for flow dosing; gas inlet

NIST integrated into rail clamp

Gas type: Oxygen O2

Inlet: NIST according to DIN EN ISO 18082

Outlet: 9/16-18 UNF

Flow Rate: 0...15 l/min

Properties

Gas Type

O2-Oxygen