

Falke with Flow Stand. O2 AGA Tubing Nipple +1AGA



Manufacturer: Greggersen

Product 910234

Number:

Description

Falke with Flow, Standard
with 1 quick-release coupling,
compact pressure regulator with cylinder connection
according to SS 367615, contains pressure gauge, incrementally adjustable
flow rate. AGA quick-release coupling at a pressure of 450 kPa.

Gas type: Oxygen / O2

Connection: Cylinder connection according to SS 367615; Connection: E1

Outlet: Tubing Nipple / 1 medium-pressure outlet type SS 87 524 30 (AGA)

Flow rate: 0...15 l/min Tubing Nipple - 120 l/min quick-release coupling

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type

O2-Oxygen