

Kolibri 2 Probe MEDAP Neo



Manufacturer: Greggersen

Product 906453

Number:

Description

KOLIBRI 2, O₂, 0...1 l/min

- Neonatal MEDAP PLUG-IN UNIT -

Flow Settings:

0 / 0.1 / 0.15 / 0.2 / 0.25 / 0.3 / 0.4 / 0.5 / 0.6 / 0.8 / 1.0 l/min

Compact Flowmeter KOLIBRI 2

with fixed orifices and selection device

Medical compact flowmeter for coupling
into a terminal unit, with handwheel - incrementally
adjustable flow rate, outlet 9/16-18 UNF thread

Ergonomically optimised adjustment device and readability
of the set flow setting possible from various viewing directions
and angles.

Frontal operating concept -> greater safety when
setting the correct flow rate

Gas type: Oxygen / O₂

Connection: MEDAP Probe

Outlet: 9/16-18 UNF

Flow Rate: 0...1 l/min

Accessories required for operation:

900619 - Hose Connection 9/16-18 UNF

Properties

Gas Type	O ₂ -Oxygen
Inlet Standard	MEDAP
Variant	Neo